



## REFINING GAS ANALYSIS

**Архангельск** (8182)63-90-72  
**Астана** (7172)727-132  
**Астрахань** (8512)99-46-04  
**Барнаул** (3852)73-04-60  
**Белгород** (4722)40-23-64  
**Брянск** (4832)59-03-52  
**Владивосток** (423)249-28-31  
**Волгоград** (844)278-03-48  
**Вологда** (8172)26-41-59  
**Воронеж** (473)204-51-73  
**Екатеринбург** (343)384-55-89  
**Иваново** (4932)77-34-06  
**Ижевск** (3412)26-03-58  
**Казань** (843)206-01-48

**Калининград** (4012)72-03-81  
**Калуга** (4842)92-23-67  
**Кемерово** (3842)65-04-62  
**Киров** (8332)68-02-04  
**Краснодар** (861)203-40-90  
**Красноярск** (391)204-63-61  
**Курск** (4712)77-13-04  
**Липецк** (4742)52-20-81  
**Магнитогорск** (3519)55-03-13  
**Москва** (495)268-04-70  
**Мурманск** (8152)59-64-93  
**Набережные Челны** (8552)20-53-41  
**Нижний Новгород** (831)429-08-12  
**Новокузнецк** (3843)20-46-81

**Новосибирск** (383)227-86-73  
**Омск** (3812)21-46-40  
**Орел** (4862)44-53-42  
**Оренбург** (3532)37-68-04  
**Пенза** (8412)22-31-16  
**Пермь** (342)205-81-47  
**Ростов-на-Дону** (863)308-18-15  
**Рязань** (4912)46-61-64  
**Самара** (846)206-03-16  
**Санкт-Петербург** (812)309-46-40  
**Саратов** (845)249-38-78  
**Севастополь** (8692)22-31-93  
**Симферополь** (3652)67-13-56  
**Смоленск** (4812)29-41-54

**Сочи** (862)225-72-31  
**Ставрополь** (8652)20-65-13  
**Сургут** (3462)77-98-35  
**Тверь** (4822)63-31-35  
**Томск** (3822)98-41-53  
**Тула** (4872)74-02-29  
**Тюмень** (3452)66-21-18  
**Ульяновск** (8422)24-23-59  
**Уфа** (347)229-48-12  
**Хабаровск** (4212)92-98-04  
**Челябинск** (351)202-03-61  
**Череповец** (8202)49-02-64  
**Ярославль** (4852)69-52-93

# REFINING GAS ANALYSIS



*Time is Money.*

Benjamin Franklin



Time is always a key factor in today's laboratories productivity.

## Master your Time with the DANI Gas Analyzers.

The ability to provide the proper configuration to meet the most challenging analytical demands comes from a long and proven experience and a deep industry knowledge. As requirements are constantly changing, even a highly reliable instrumentation could not be enough to succeed in getting trustworthy results: complete and guaranteed solutions are essential to comply with the latest industry standards and specifications.

After a long working relationship with its customers to know and to best match their real needs, DANI Instruments has developed key analytical solutions that cover a broad array of applications, requirements and protocols in the environmental industry.

Master DANI Air Analyzers are **PRE-CONFIGURED, PRE-ASSEMBLED AND FACTORY-TESTED SYSTEMS** specifically designed for peculiar analyses. The analyzers include the **HARDWARE**, the **SOFTWARE**, **COLUMNS AND CONSUMABLES**, the **OPTIMIZED ANALYSIS METHOD**, the **ANALYTICAL CONDITIONS**, and the **DOCUMENTATION** to run up your analysis from day one.

### **PRE-CONFIGURED, PRE-ASSEMBLED AND FACTORY-TESTED SYSTEMS**

The installation process is faster than ever before and all the startup procedure is oversimplified ensuring immediate analytical performance and results.

### **HARDWARE AND SOFTWARE**

DAA Analyzers are pre-engineered systems based on the versatility, flexibility and robustness of the proven Master GC hardware. All the Master GC parameters are set prior the shipment.

### **COLUMNS AND CONSUMABLES**

No more doubts about the proper column, parts and supplies. DAA Analyzers are delivered with all you may need for your analysis.\*

### **OPTIMIZED ANALYSIS METHOD**

Analytical methods are pre-loaded to be immediately used for the determination of pollutants in air. Whenever possible, reduction of analysis time and amounts of toxic solvents are considered. Method development time and costs are thus dramatically reduced.

### **ANALYTICAL CONDITIONS**

DAA Analyzers are designed to perfectly accomplish the analytical conditions of interest.

### **DOCUMENTATION**

A getting started manual, calibration and method files, and all the information for a quick startup are included.



\* Chemicals are not supplied




## DANI MASTER GC

- Intuitive and easy-to-use touchscreen interface.
- Complete range of injectors:
  - Split/Splitless Injector
  - Programmable Temperature Vaporizer (PTV)
  - Packed Injector
- Wide selection of detectors:
  - Flame Ionization Detector (FID)
  - Electron Capture Detector (ECD)
  - Nitrogen-Phosphorus Detector (NPD)
  - Flame Photometric Detector (FPD)
  - Thermal Conductivity Detector (TCD)
  - Pulse Discharge Detector (PDD)
  - Master TOF-MS Time of Flight Mass Spectrometer
- Extensive Choice of Dedicated Devices:
  - Auxiliary ovens
  - Gas sampling and switching valves
  - Liquid sampling valves
  - GC oven cryogenic cooling device
  - Methanizer

## DANI MASTER AUX

### Auxiliary Ovens for Isothermal Temperature

- Up to 7 valves (5 gas + 2 liquids)
- Up to 250° C
- Up to 2 auxiliary gas valves
- Up to 6 needle valves
- Extremely compact design
- Two models available with different capacities to house valves and columns for dedicated analyses



## DANI MASTER AS


### Liquid Autosampler

- Superior flexibility, repeatability, and performance
- Easy-to-use
- 160 samples capacity
- No sample degradation or solvent evaporation
- up to 7 syringe capacity types
- unmatched injection capabilities

## DANI MASTER SHS

### Static Headspace Sampler

- A robust and flexible system to meet complex and versatile needs
- Reliable results and exceptional reproducibility
- Intuitive, powerful, and straightforward user interface
- Highest sample capacity with a 120-position vial tray
- Unlimited priority sample position
- Unmatched oven capacity: 18 vials simultaneously



# REFINING GAS ANALYSIS

## TABLE OF CONTENTS

Configuration	Typical Detected Compounds	Official Methods Met	Sample Type
<b>DGA-003 Refining Gas Analyzer for Permanent Gas, H<sub>2</sub>, and Hydrocarbons</b>			
Three channels with FID, dual $\mu$ TCD, micro-packed columns, three valves, and auxiliary oven.	hydrogen helium nitrogen carbon dioxide methane ethane	propane butane <i>iso</i> -butane pentane C <sub>6</sub> + (hexane and higher hydrocarbons)	Gas Liquid
		ASTM 1945 ASTM 1946 GPA 2177 ASTM 2163 UOP 539 ISO 6874	
<b>DGA-010 Refining Gas Analyzer for Hydrogen and Helium</b>			
Dual channels with $\mu$ TCD, micro-packed columns, two valves, auxiliary oven.	hydrogen helium		Gas
		ASTM 2504	
<b>DGA-014 Refining Gas Analyzer for Natural Gas with Permanent Gas separation - H<sub>2</sub> and He</b>			
Dual channels with dual $\mu$ TCD, three valves, auxiliary oven.			
		ASTM 1945 ASTM 1946 GPA 2177 ASTM 2504	
<b>DGA-021 Refining Gas Analyzer for Gas and Liquefied Hydrocarbons</b>			
Dual channels with FID, two gas valves and one liquid valve, capillary columns, auxiliary oven.	methane ethane propane <i>iso</i> -butane	<i>n</i> -butane <i>iso</i> -pentane <i>n</i> -pentane hexane	Gas Liquid
		ASTM 2163	
<b>DGA-028 Refining Gas Analyzer for Permanent Gas, Hydrocarbons, H<sub>2</sub>, and CO<sub>2</sub></b>			
Three channels with $\mu$ TCD, FID, three valves, auxiliary oven.	nitrogen carbon dioxide hydrogen methane ethane	propane <i>i</i> -butane <i>n</i> -butane <i>i</i> -pentane <i>n</i> -pentane <i>n</i> -hexane	Gas
		ASTM 1945 ASTM 1946 GPA 2177	
<b>DGA-033 Refining Gas Analyzer for Permanent Gas, CO<sub>2</sub>, Hydrocarbons, COS, H<sub>2</sub>S, CS<sub>2</sub></b>			
Dual channels with FID, $\mu$ TCD, FPD, three valves, and auxiliary oven.	oxygen nitrogen carbon dioxide carbonyl sulfide hydrogen sulfide carbon disulfide methane	ethane propane <i>i</i> -butane <i>n</i> -butane <i>i</i> -pentane <i>n</i> -pentane <i>n</i> -hexane	Gas
<b>DGA-034 Refining Gas Analyzer for Permanent Gas and Hydrocarbons</b>			
Two channels with TCD, packed columns, two valves, and auxiliary oven.	oxygen nitrogen methane carbon dioxide ethane	propane butane hexane	Gas
		ASTM 1945 ASTM 1946 GPA 2177	
<b>DGA-035 Refining Gas Analyzer for Permanent Gas, H<sub>2</sub>, and Hydrocarbons</b>			
Two channels with FID, $\mu$ TCD, micro-packed and capillary columns, two valves for gas and one valve for liquids, auxiliary oven.	hydrogen oxygen nitrogen carbon dioxide methane carbon monoxide	ethane propane butane pentane C <sub>6</sub> + (hexane and higher hydrocarbons)	Gas Liquid
		ASTM 1945 ASTM 1946 GPA 2177	

Configuration	Typical Detected Compounds	Official Methods Met	Sample Type
<b>DGA-036 Refining Gas Analyzer for Permanent Gas, H<sub>2</sub>, and Hydrocarbons</b> Three channels with $\mu$ TCD, micro-packed columns, three valves, auxiliary oven.	hydrogen oxygen nitrogen carbon dioxide methane carbon monoxide	ethane propane butane pentane C <sub>6</sub> + (hexane and higher hydrocarbons)	ASTM 1945 ASTM 1946 GPA 2177  Gas

# REFINING GAS ANALYSIS

## DGA-003

### Natural Gas Analyzer for Permanent Gas, H<sub>2</sub> and Hydrocarbons Determination

DGA-003 is a three-channels system that allows the accurate determination of permanent gases, H<sub>2</sub>, and hydrocarbons through the two channels equipped with the  $\mu$ TCD detectors and the channel equipped with the FID detector. The proper separation of permanent gases and hydrocarbons is possible thanks to the presence of a micro-packed Hayesep column for each TCD channel.

DANI DGA-003 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, H<sub>2</sub> and Hydrocarbons.

**Master DGA-003**

#### SAMPLE TYPE:

Gas    Liquid

#### TYPICAL DETECTED COMPOUNDS:

H<sub>2</sub>, He, O<sub>2</sub>, N<sub>2</sub>, CO<sub>2</sub>, CH<sub>4</sub>, ethane, propane, butane, *iso*-butane, pentane, and C<sub>6</sub>+ (hexane and higher hydrocarbons)

#### OFFICIAL METHODS MET:

ASTM 1945, ASTM 1946, GPA2177, ASTM 2163, UOP 539

#### KEY FEATURES & BENEFITS:


5 micro-packed columns installed for a complete and accurate separation;

Gas sample is introduced simultaneously into each analytical channel by the 10-port valve;


The external valve gives the possibility to introduce liquid samples into the third channel;

The TCD detector has advantages for speed, accuracy and reduced dead volumes for permanent gases and H<sub>2</sub>;


The FID detector is very sensitive to the hydrocarbons components with an excellent linearity across a wide range of concentrations;



µTCD chromatogram of the gas standard mixture. It is possible to recognize the peaks of CO<sub>2</sub>, C<sub>2</sub>H<sub>6</sub>, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, C<sub>3</sub>H<sub>8</sub>



µTCD (channel B) chromatogram of the gas standard mixture (C5). Labeled peak correspond to: H<sub>2</sub>. All the other compounds of the standard mix are back flushed through Hayasep Q since they are not of interest for this channel and to avoid that they enter in the Molsieve column.



FID chromatogram of the gas standard mixture (C6). Labeled peaks correspond to: CH<sub>4</sub>, ethane, propane, isobutane, n-butane, isopentane, n-pentane, and n-hexane.

EXPERIMENTAL PARAMETERS		
<b>Master GC Analyzer</b>		
Columns	Channel 1	Hayesep Q (2m,1/16",80/100 mesh) Molsieve 5A (2m,1/16",80/100 mesh)
	Channel 2	Hayesep Q (2m,1/16",80/100 mesh) Molsieve 5A (2m,1/16",80/100 mesh)
	Channel 3	Alumina (1m,1/16",80/100 mesh)
<b>GC Oven</b>		
Temperature	32°C (for 10 min), 50°C/min to 150°C (for 20 min)	
<b>Injector A : PK</b>		
Temperature	40°C	
Carrier Gas	Helium	
Maximum Pressure	n/a	
Flow	6 ml/min	
<b>Injector B : PK</b>		
Temperature	40°C	
Carrier Gas	Nitrogen	
Maximum Pressure	n/a	
Flow	8 ml/min	
<b>Injector C: SL/IN</b>		
Temperature	250°C	
Carrier Gas	Helium	
Split Ratio	1:10	
Flow	12 ml/min	
<b>Detector A : TCD Channel 1</b>		
Temperature	150°C	
Main filament temp.	220°C	
Main filament safety	Injector A	
Min Half-Peak Width	1s	
Digital Acquisition Rate	10 Hz	
Signal zeroing	50 mV	
<b>Detector B : TCD Channel 2</b>		
Temperature	150°C	
Main Filament Temp.	220°C	
Main Filament Safety	Injector B	
Min. Half-Peak Width	1s	
Digital Acquisition Rate	10 Hz	
Signal Zeroing	50 mV	
<b>Detector C: FID</b>		
Temperature	150°C	
AUX Flow	25 ml/min	
H <sub>2</sub> Flow	40 ml/min	
Air Flow	280 ml/min	
<b>Auxiliary Pressure</b>		
AUX 1 Pressure	1.5 bar	
AUX 2 Pressure	2.2 bar	
<b>Auxiliary Temperature</b>		
Aux temp - Aux oven temp	80°C	

# REFINING GAS ANALYSIS

## DGA-010

### Natural Gas Analyzer for Hydrogen and Helium

The DGA-010 Analyzer enables the analysis of hydrogen and helium in a gas mixture which is performed using a Packed (PK) injector, micro-packed columns and a  $\mu$ TCD detector.

The instrument is equipped with two identical channels with this same configuration.

Moreover, the gas chromatograph is configured with two 8-port valves located in the GC auxiliary oven and a backflush system.

DANI DGA-010 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Hydrogen and Helium

*Master* **DGA-010**

#### SAMPLE TYPE:

Gas

#### TYPICAL DETECTED COMPOUNDS:

H<sub>2</sub> and He

#### OFFICIAL METHODS MET:

ASTM 2504


#### KEY FEATURES & BENEFITS:

With the DGA-010 Analyzer it is possible to determine Hydrogen and Helium in a hydrocarbon mixture and separate them from the hydrocarbons using a backflush;

Valve system configuration and their default positions of switching are factory set up;

All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;





µTCD channel A chromatogram of the Gas standard mixture. It is possible to see the peaks of Helium and H2.

EXPERIMENTAL PARAMETERS FOR H <sub>2</sub> AND HE ANALYSIS CHANNEL A	
<b>Master GC Analyzer</b>	
Columns	Channel A Hayesep Q (1m, 1mm ID, 1/16", 100/120 mesh) Molsieve (2m, 1mm ID, 1/16", 80/100 mesh)
<b>GC Oven</b>	
Temperature	at 3 min -33°C to 60°C at 20°C/min for 5 min
<b>Injector A: PK</b>	
Temperature	200°C
Carrier Gas	Helium
Maximum Pressure	n/a
Flow	15 ml/min
<b>Detector A : µTCD</b>	
Temperature	150°C
Filament Temperature	210°C
Signal Zeroing	50 mV
Filament Safety	Injector A
<b>Auxiliary Pressure</b>	
AUX 1 Pressure	2.7 bar
<b>Auxiliary Temperature</b>	
Aux temp 1 - Aux oven temp	80°C



µTCD channel B chromatogram of the Gas standard mixture. It is possible to see the peaks of Helium and H2.

EXPERIMENTAL PARAMETERS FOR H <sub>2</sub> AND HE ANALYSIS CHANNEL B	
<b>Master GC Analyzer</b>	
Columns	Channel B Hayesep Q (1m, 1mm ID, 1/16", 100/120 mesh) Molsieve (2m, 1mm ID, 1/16", 80/100 mesh)
<b>GC Oven</b>	
Temperature	at 3 min -33°C to 60°C at 20°C/min for 5 min
<b>Injector B: PK</b>	
Temperature	200°C
Carrier Gas	Helium
Maximum Pressure	n/a
Flow	15 ml/min
<b>Detector B : µTCD</b>	
Temperature	150°C
Filament Temperature	210°C
Signal Zeroing	50 mV
Filament Safety	Injector B
<b>Auxiliary Pressure</b>	
AUX 1 Pressure	4 bar
<b>Auxiliary Temperature</b>	
Aux temp 1 - Aux oven temp	80°C

# REFINING GAS ANALYSIS

## DGA-014

### Natural Gas Analyzer for NGA with Permanent Gas Separation - H<sub>2</sub> and He

The the DGA-014 Analyzer is a dual-channel system for the analysis of Permanent Gases, Hydrogen and Helium. The first one, for permanent gas analysis, is configured with a packed injector (PK), micro-packed columns and a Micro Thermal Conductivity Detector ( $\mu$ TCD).

The analysis of Hydrogen in a gas mixture, on the second channel, is performed by a Packed Injector (PK), micro-packed column and the Micro Thermal Conductivity Detector ( $\mu$ TCD).

DANI DGA-014 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, H<sub>2</sub> and He.

*Master* **DGA-014**

#### SAMPLE TYPE:

#### TYPICAL DETECTED COMPOUNDS:

#### OFFICIAL METHODS MET:

ASTM 1945, ASTM 1946, GPA 2177, ASTM 2504


#### KEY FEATURES & BENEFITS:

With the DGA-014 it is possible to determine O<sub>2</sub>, CH<sub>4</sub>, N<sub>2</sub>, CO<sub>2</sub>, C<sub>2</sub> and higher hydrocarbons in Natural Gas Samples and separate the components above C<sub>3</sub> using a backflush;

The use of an optimized combination of columns allows the excellent separation of all the components and more reproducible characteristics;


The  $\mu$ TCD detector guarantees linear dynamic range, speed, reduced dead volumes and precision;

All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;



uTCD A chromatogram of the Hydrocarbons standard mixture. It is possible to recognize the peaks of CO<sub>2</sub>, C<sub>2</sub>H<sub>6</sub>, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, C<sub>3</sub>H<sub>8</sub>

EXPERIMENTAL PARAMETERS FOR NATURAL GAS ANALYSIS CHANNEL A		
<b>Master GC Analyzer</b>		
Columns	Channel A	Hayesep T (1m, 1/16", 80/100 mesh) Hayesep Q (1m, 1/16", 100/120 mesh) Molsieve 5A (2m, 1/16, 80/100 mesh)
<b>GC Oven</b>		
Temperature	at 15 min 33°C to 150°C at 20°C/min for 10 min	
<b>Injector A: PK</b>		
Temperature	200°C	
Carrier Gas	Helium	
Maximum Pressure	n/a	
Flow	15 ml/min	
<b>Detector A : μTCD</b>		
Temperature	150°C	
Filament Temp	210°C	
Signal Zeroing	50 mV	
Filament Safety	Injector A	
<b>Auxiliary Temperature</b>		
Aux Temp 1 - Aux Oven Temp	80°C	
<b>Auxiliary Pressure</b>		
AUX Gas 1	3.2 bar	



uTCD channel B chromatogram of the Hydrogen and Hydrocarbons standard mixture. It is possible to see the peaks of Helium and H<sub>2</sub>

EXPERIMENTAL PARAMETERS FOR H <sub>2</sub> AND HE ANALYSIS CHANNEL B		
<b>Master GC Analyzer</b>		
Columns	Channel B	Hayesep Q (1m, 1mm ID, 1/16", 100/120 mesh) Molsieve (2m, 1mm ID 1/16", 80/100 mesh)
<b>GC Oven</b>		
Temperature	at 4 min 33°C to 150°C at 20°C/min for 10 min	
<b>Injector B: PK</b>		
Temperature	200°C	
Carrier Gas	Helium	
Maximum Pressure	n/a	
Flow	17.2 ml/min	
<b>Detector B : μTCD</b>		
Temperature	150°C	
Filament Temperature	210°C	
Signal Zeroing	50mV	
Filament Safety	Injector B	
<b>Auxiliary Temperature</b>		
Aux temp 1 - Aux oven temp	80°C	
<b>Auxiliary Pressure</b>		
AUX 3 Pressure	3 bar	

# REFINING GAS ANALYSIS

## DGA-021

### Refining Gas Analyzer for Gas and Liquefied Hydrocarbons

The analysis of hydrocarbons in a gas or liquid phase is performed using a split/splitless injector, capillary column and a Flame Ionization Detector.

The GC is configured with two valves for gas sampling and equilibrating with atmosphere pressure: two 6-port valves are located in the GC auxiliary oven. Moreover, a 4-port valve for liquid sampling is available.

DANI DGA-021 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Gas and Liquefied Hydrocarbons.

*Master* **DGA-021**

#### SAMPLE TYPE:

Gas    Liquid

#### TYPICAL DETECTED COMPOUNDS:

CH<sub>4</sub>, ethane, propane, iso- and n-butane, iso- and n-pentane, hexane;


#### OFFICIAL METHODS MET:

ASTM 2163

#### KEY FEATURES & BENEFITS:


Hydrocarbons analysis is performed using the two 6-port valves system, a capillary column and an FID detector;

All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;



FID chromatogram of the gas standard mixture. It is possible to recognize the peaks of CH<sub>4</sub>, C<sub>2</sub>H<sub>6</sub>, C<sub>3</sub>H<sub>8</sub>, i-butane, n-butane, i-pentane, n-pentane

EXPERIMENTAL PARAMETERS FOR HYDROCARBONS DETERMINATION	
<b>Master GC Analyzer</b>	
Columns	Al <sub>2</sub> O <sub>3</sub> /Na <sub>2</sub> SO <sub>4</sub> (50 m, 0.25 mm ID, 4µm)
<b>GC Oven</b>	
Temperature	At 0.00 min from 35°C to 180°C at 20°C/min - for 7 min
<b>Injector B: SL/IN</b>	
Temperature	20°C
Carrier Gas	Helium
Split Ratio	1:50
Flow	1 ml/min
<b>Detector B: FID</b>	
Temperature	300°C
Aux Flow (N <sub>2</sub> )	25 ml/min
H <sub>2</sub> Flow	40 ml/min
Air Flow	280 ml/min
<b>Aux Temperature</b>	
Aux Temp 1 (°C) - Aux Oven Temp	80°C



FID chromatogram of the liquid hydrocarbons mixture. A real sample of lighter fuel was analyzed. It was possible to identify the peaks of CH<sub>4</sub>, C<sub>2</sub>H<sub>6</sub>, C<sub>3</sub>H<sub>8</sub>, iso and n-butane, iso and n-pentane.

EXPERIMENTAL PARAMETERS FOR LIQUID HYDROCARBONS DETERMINATION	
<b>Master GC Analyzer</b>	
Columns	Al <sub>2</sub> O <sub>3</sub> /Na <sub>2</sub> SO <sub>4</sub> (50 m, 0.25 mm ID, 4µm)
<b>GC Oven</b>	
Temperature	At 0.00 min from 35°C to 200°C at 20°C/min - for 7 min
<b>Injector B: SL/IN</b>	
Temperature	20°C
Carrier Gas	Helium
Split Ratio	1:20
Flow	1 ml/min
<b>Detector B: FID</b>	
Temperature	300°C
Aux Flow (N <sub>2</sub> )	25 ml/min
H <sub>2</sub> Flow	40 ml/min
Air Flow	280 ml/min
<b>Aux Temperature</b>	
Aux Temp 1 (°C) - Aux Oven Temp	80°C

# REFINING GAS ANALYSIS

## DGA-028

### Refining Gas Analyzer for Permanent Gas, Hydrocarbons, H<sub>2</sub> and CO<sub>2</sub>

The DGS-028 Analyzer enables the analysis of permanent gas, CO<sub>2</sub> and hydrocarbons up to C<sub>6</sub> using three different pathways, channel 1, channel 2 and channel 3 equipped as follows:

Channel 1 is composed by a flow control module injector, a 10-port valve for sampling and switching, a 6-port valve for switching, two types of columns: Plot Q and Molsieve 5A and a thermal conductivity detector (TCD). This channel fullfills the need of an accurate determination of permanent gas, CO<sub>2</sub> and CH<sub>4</sub>.

Channel 2 is composed by a flow control module injector, a 10-port valve for sampling and switching, two types of columns: Plot Q and Molsieve 5A and a thermal conductivity detector (TCD). This channel fullfills the need of an accurate determination of H<sub>2</sub>.

Channel 3 is composed by a SL/IN injector, a 6-port valve for switching, an Alumina column and a flame ionized detector (FID). Sampling is made by a single "sample in" for all channels.

DANI DGA-028 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, Hydrocarbons, H<sub>2</sub> and CO<sub>2</sub>

**Master DGA-028**

#### SAMPLE TYPE:

Gas

#### TYPICAL DETECTED COMPOUNDS:

O<sub>2</sub>, N<sub>2</sub>, CO<sub>2</sub>, H<sub>2</sub>, methane, ethane, propane, *i*-butane, *n*-butane, *i*-pentane, *n*-pentane, *n*-hexane

#### OFFICIAL METHODS MET:


ASTM 1945, ASTM 1946, GPA 2177

#### KEY FEATURES & BENEFITS:


Three channel system for a comprehensive, ast, and accurate analysis of permanent gases, hydrocarbons, CO<sub>2</sub>, H<sub>2</sub>, and CH<sub>4</sub>;

Higher sensitivity and precision are obtained through the use of the μTCD and FID detectors;


All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;



µTCD chromatogram of gas standard mixture of permanent gases and CH4



µTCD chromatogram of gas standard mixture



FID chromatogram of the gas standard mixture

**EXPERIMENTAL PARAMETERS FOR PERMANENT GAS, CO<sub>2</sub>, H<sub>2</sub> AND HYDROCARBONS**

Master GC Analyzer		
Columns	Channel 1	Porapak Q (3m, 80/100 mesh, 1/8") Molsieve 5A (3m, 80/100 mesh, 1/8")
	Channel 2	Porapak Q (3m, 80/100 mesh, 1/8") Molsieve 5A (3m, 80/100 mesh, 1/8")
	Channel 3	Alumina (1m, 80/100 mesh, 1/16")
GC Oven		
Temperature	1) T : 40°C, Time: 15 min, Rate: 10°C/min 2) T: 200°C, Time: 4 min	
Injector A: SL/IN		
Temperature	250°C	
Carrier Gas	Helium	
Split Ratio	1:80	
Flow	1 ml/min	
Injector B: SL/IN		
Carrier Gas	Helium	
Flow	8 ml/min	
Injector C: SL/IN		
Carrier Gas	Helium	
Flow	20 ml/min	
Det A: TCD		
Temperature Control	250°C	
Filaments Control	7.40 V	
Polarity	Negative	
Max Current	180 mA	
Filament Safety	Aux Gas 3	
Filament Safety	Aux Gas 1	
Min. Half-Peak Width	0.60 s	
Digital Acquisition Rate	25 Hz	
Signal Zeroing	10 mV	
Det B: TCD		
Temperature Control	250°C	
Filaments Control	7.40 V	
Polarity	Positive	
Max Current	180 mA	
Filament Safety	Inj C	
Filament Safety	Aux Gas 2	
Min. Half-Peak Width	0.60 s	
Digital Acquisition Rate	25 Hz	
Signal Zeroing	10 mV	
Det C: FID		
Temperature	300°C	
AUX flow	25 ml/min	
H <sub>2</sub> flow	40 ml/min	
Air flow	280 ml/min	
Auxiliary Pressure		
AUX 1 pressure	0.70 bar	
AUX 2 pressure	1.10 bar	
AUX 3 pressure	3.00 bar	
Aux Temperature		
Aux Temperature	80°C	

# REFINING GAS ANALYSIS

## DGA-033

Refinery Gas Analyzer for Permanent Gas, CO<sub>2</sub>, Hydrocarbons, COS, H<sub>2</sub>S and CS<sub>2</sub>

This MASTER GC Analyzer enables the analysis of permanent gas, CO<sub>2</sub>, COS, H<sub>2</sub>S, CS<sub>2</sub> and hydrocarbons up to C<sub>6</sub>.

The analysis are performed with two different pathway equipped as follows:

Channel 1 is composed by a PK injector, two types of columns (Plot Q and Molesieve), a micro Thermal Conductivity Detector ( $\mu$ TCD), a Flame Photometric Detector (FPD);

Channel 2 is composed by a SL/IN injector, a Plot Q column and a Flame Ionized Detector (FID).

Moreover, the gas chromatograph is configured with a 10-ports valve for sampling, two 6-ports valves for switching and a 8-ports liquid valve all located in the GC auxiliary oven.

DANI DGA-033 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, CO<sub>2</sub>, Hydrocarbons, COS, H<sub>2</sub>S, and CS<sub>2</sub>

**Master DGA-033**

### SAMPLE TYPE:

Gas

### TYPICAL DETECTED COMPOUNDS:

O<sub>2</sub>, N<sub>2</sub>, CO<sub>2</sub>, COS, H<sub>2</sub>S, CS<sub>2</sub>, methane, ethane, propane, *i*-butane, *n*-butane, *i*-pentane, *n*-pentane, *n*-hexane;

### OFFICIAL METHODS MET:

### KEY FEATURES & BENEFITS:


Typical components of interest in this type of gas analysis are: O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, CO, COS, H<sub>2</sub>S, CS<sub>2</sub>, methane, ethane, propane, *i*-butane, *n*-butane, *i*-pentane, *n*-pentane and *n*-hexane.

With this analyzer it is possible to separate and analyze these compounds or part of them.


Higher sensitivity and precision are obtained through the use of the  $\mu$ TCD and FID detectors;

All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;






µTCD (Chn 1) chromatogram of the gas standard mixture



FPD (Chn 1) analysis of sulphur compounds mixture



FID (Chn 2) analysis of C6 hydrocarbons mixture

**EXPERIMENTAL PARAMETERS FOR PERMANENT GAS, CO<sub>2</sub>, HYDROCARBONS AND SULPHUR COMPOUNDS**

Master GC Analyzer		
Columns	Channel 1	Plot Q (3m, 1/16 OD) Molsieve (3m, 1/16 OD)
	Channel 2	Plot Q (3m, 1/16 OD)
GC Oven		
Temperature	1) T : 35°C, Time: 15 min, Rate: 30°C/min 2) T: 200°C, Time: 11 min	
Injector A: PK		
Temperature	250°C	
Carrier Gas	Helium	
Flow	10 ml/min	
Injector B: SL/IN		
Temperature	250°C	
Carrier Gas	Helium	
Split Ratio	1:25	
Flow	8 ml/min	
Det A: FPD		
Temperature Control DBB	250°C	
Detector Head Temperature (Aux Temp 2)	130°C	
H <sub>2</sub> Flow Rate	200 ml/min	
Air 1 Flow Rate	160 ml/min	
Air 2 Flow Rate	160 ml/min	
Aux Flow Rate (N <sub>2</sub> )	15 ml/min	
Photomultiplier voltage	0.72 kV	
Range	1	
Min. Half-Peak Width	0.60 s	
Digital Acquisition Rate	25 Hz	
Det B: FID		
Temperature	300°C	
Aux Flow Rate (N <sub>2</sub> )	25 ml/min	
H <sub>2</sub> Flow Rate	40 ml/min	
Air Flow Rate	280 ml/min	
Range	10	
Min. Half Peak Width	0.60 s	
Digital Acq. Rate	25 Hz	
Det C: µTCD		
Temperature Control	160°C	
Main Filament Temperature	190°C	
Main Filament Safety	Injector A	
Min. Half-Peak Width	0.60 s	
Digital Acquisition Rate	25 Hz	
Signal Zeroing	10 mV	
Auxiliaries		
Aux Temp 1- for Aux Oven	80°C	
Aux Temp 2 - For Detector Head temp FPD	130°C	
Aux Gas (N <sub>2</sub> ) - for Backflush	1.7 bar	

# REFINING GAS ANALYSIS

## DGA-034

### Refinery Gas Analyzer for Permanent Gas and Hydrocarbons

The analysis with the DGA-034 Analyzer is performed using two Packed Column Injectors (PK), three micropacked columns and two Thermal Conductivity Detectors (TCD).

Moreover, the Gas Chromatograph is configured with two valves located in the auxiliary oven: a 10-port valve for injection in the analytical channels and a 6-port valve on Channel 1.

DANI DGA-034 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, H<sub>2</sub> and Hydrocarbons

*Master* **DGA-034**

#### SAMPLE TYPE:

Gas

#### TYPICAL DETECTED COMPOUNDS:

O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, CO<sub>2</sub>, ethane, propane, butane, pentane, and hesane;

#### OFFICIAL METHODS MET:

ASTM 1945, ASTM 1946, GPA 2177;


#### KEY FEATURES & BENEFITS:

Gas analysis is performed using this complete system that allows to determinate permanent gas and hydrocarbons using the channel (channel 1) equipped with TCD, a packed Plot Q column and a packed Molsieve column, a 10-port valve and a 6-port valve;


For the determination of hydrocarbons, Channel 2 is equipped with the TCD, the same 10-port valve for injection and a packed Plot Q column;

The gas sample is introduced simultaneously to each analytical channel through the 10-port valve equipped with sampling loops:

All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;



Chromatogram of the gas standard mixture on Channel 1.  
The peaks correspond to: CO<sub>2</sub>, C<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, C<sub>3</sub>, iso-C<sub>4</sub>, n-C<sub>4</sub>, iso-C<sub>5</sub>, n-C<sub>5</sub>, C<sub>6</sub>



Chromatogram of the gas standard mixture on Channel 2.  
The peaks correspond to: O<sub>2</sub> + N<sub>2</sub>, CH<sub>4</sub>, CO<sub>2</sub>, C<sub>2</sub>, C<sub>3</sub>, iso-C<sub>4</sub>, n-C<sub>4</sub>, iso-C<sub>5</sub>, n-C<sub>5</sub>, C<sub>6</sub>

EXPERIMENTAL PARAMETERS FOR PERMANENT GAS, CO <sub>2</sub> , HYDROCARBONS AND SULPHUR COMPOUNDS		
<b>Master GC Analyzer</b>		
Columns	Channel 1	<b>Hayesep Q</b> (2,5m, 1/8", 80/100 mesh) <b>Molsieve 5A</b> (2.5m, 1/8", 80/100mesh)
	Channel 2	<b>Hayesep Q</b> (2m, 1/16", 80/100mesh)
<b>GC Oven</b>		
Temperature	1) T: 40°C, Time: 20 min, Rate: 20°C/min 2) T: 200°C, Time: 8 min	
<b>Injector A: PK</b>		
Temperature	250°C	
Carrier Gas	Helium	
Column Type	Packed - Max Pressure: 3 bar	
Flow	10 ml/min	
<b>Injector B: PK</b>		
Temperature	250°C	
Carrier Gas	Helium	
Column Type	Packed - Max pressure: 2 bar	
Flow	12 ml/min	
<b>Channel 1 - Det A: TCD</b>		
Temperature	250°C	
Voltage	7.20 V	
Polarity	Negative	
Max Current	180 mA	
Filament Safety	Inj A	
Filament Safety	Aux Gas 1	
Range	1	
Min. Half-Peak Width	0.60 s	
Digital Acq. Rate	25 Hz	
Signal Target	10 mV	
<b>Channel 1 - Det B: TCD</b>		
Temperature	250°C	
Voltage	7.30 V	
Polarity	Negative	
Max Current	180 mA	
Filament Safety	Inj B	
Filament Safety	Aux Gas 2	
Range	1	
Min. Half-Peak Width	0.60 s	
Digital Acq. Rate	25 Hz	
Signal Target	10 mV	
<b>Auxiliary Temperature</b>		
Aux Temp (°C) - Aux Oven Temp	120°C	
<b>Auxiliary Pressure</b>		
Aux Press 1	0.65 bar	
Aux Press 2	0.50 bar	

# REFINING GAS ANALYSIS

## DGA-035

### Refinery Gas Analyzer for Permanent Gas, H<sub>2</sub> and Hydrocarbons

The analysis with the DGA-035 Analyzer is performed using two injector flow controls, a split/splitless injector, micro-packed columns, a micro Thermal Conductivity Detector ( $\mu$ TCD) and a Flame Ionization Detector (FID). Moreover, the gas chromatograph is configured with three valves: an external 4-ports valve (for introduce liquid samples in the system) plus a 10-ports valve and a 6-ports valve located in the auxiliary oven.

DANI DGA-035 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, H<sub>2</sub> and Hydrocarbons

*Master* **DGA-035**

#### SAMPLE TYPE:

Gas    Liquid

#### TYPICAL DETECTED COMPOUNDS:

H<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CO<sub>2</sub>, CH<sub>4</sub>, CO, ethane, propane, butane, pentane and C<sub>6</sub>+ (hexane and higher hydrocarbons)

#### OFFICIAL METHODS MET:

ASTM 1945, ASTM 1946, GPA 2177;


#### KEY FEATURES & BENEFITS:

Gas analysis is performed using this complete system that allows to determinate permanent gases, H<sub>2</sub> and hydrocarbons using the channel equipped with  $\mu$ TCD detector and the channel equipped with FID detector.


The separation of permanent gases and hydrocarbons is made possible by the presence of a micro-packed Hayesep Q column and a micro-packed Molsieve column for each channel with  $\mu$ TCD.

The gas sample is introduced simultaneously to each analytical channel (Chn 1 and 2) by the 10-ports valve equipped with two sample loops. Moreover the external valve gives the possibility to introduce liquid samples in Channel 3.


All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;



mTCD (Channel 1) chromatogram of the gas standard mixture.  
 Labelled peaks correspond to: O<sub>2</sub>, CO<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, C<sub>2</sub>, C<sub>3</sub>, iso-C<sub>4</sub>, n-C<sub>4</sub>, iso-C<sub>5</sub>, n-C<sub>5</sub>, C<sub>6</sub>



mTCD (Channel 2) chromatogram of the gas standard mixture.  
 Labelled peaks correspond to: H<sub>2</sub> and CH<sub>4</sub>



FID chromatogram of the real sample of Butane

EXPERIMENTAL PARAMETERS FOR PERMANENT GAS, CO <sub>2</sub> , HYDROCARBONS AND SULPHUR COMPOUNDS		
<b>Master GC Analyzer</b>		
Columns	Channel 1	Hayesep Q (3m, 1/16", 80/100 mesh) Molsieve 5A (3m, 1/16", 80/100mesh)
	Channel 2	Hayesep Q (3m, 1/16", 80/100mesh) Molsieve 5A (3m, 1/16", 80/100mesh)
	Channel 2	Alumina (50m, 0.32mm)
<b>GC Oven</b>		
Temperature	35°C (for 11 min) at 30°C/min to 200°C (for 13 min)	
<b>Injector A: PK</b>		
Temperature	OFF	
Carrier Gas	Helium	
Column Type	Capillary - 200m, 0.53mm	
Flow	20 ml/min	
<b>Injector B: PK</b>		
Temperature	OFF	
Carrier Gas	Nitrogen	
Column Type	Capillary . 200m, 0.53mm	
Flow	10 ml/min	
<b>Inj C : SL/IN</b>		
Temperature	200°C	
Carrier Gas	Helium	
Split Ratio	1:50	
Flow	4ml/min	
<b>Det A : μTCD Channel 1</b>		
Temperature	200°C	
Main Filament Temperature	190°C	
Filament Safety	Injector A	
Min. Half-Peak Width	1s	
Digital Acquisition Rate	10 Hz	
Signal Zeroing	10 mV	
<b>Det B : μTCD Channel 2</b>		
Temperature	160°C	
Main Filament Temperature	190°C	
Filament Safety	Injector B	
Min Half-Peak Width	1 s	
Digital Acquisition Rate	10 Hz	
Signal Zeroing	10 mV	
<b>Det C : FID</b>		
Temperature	300°C	
AUX Flow	25ml/min	
H2 Flow	40ml/min	
Air Flow	280ml/min	
<b>Auxiliary Pressure</b>		
Aux Press 1	2.6 bar	
Aux Press 2	1.45 bar	
<b>Auxiliary Temperature</b>		
Aux Temp (°C) - AUx Oven Temp	120°C	

# REFINING GAS ANALYSIS

## DGA-036

### Refinery Gas Analyzer for Permanent Gas, H<sub>2</sub> and Hydrocarbons

The analysis with MASTER GC Analyzer is performed using three injector flow controls, micro-packed columns and two micro Thermal Conductivity Detectors (mTCD).

Moreover, the gas chromatograph is configured with a solenoid valve and three valves located in the auxiliary oven: a 8-ports valve for Channel 1, plus a 10-ports valve and a 6-ports valve for Channel 2.

DANI DGA-036 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Permanent Gas, H<sub>2</sub> and Hydrocarbons

*Master* **DGA-036**

#### SAMPLE TYPE:

Gas

#### TYPICAL DETECTED COMPOUNDS:

H<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CO<sub>2</sub>, CH<sub>4</sub>, CO, ethane, propane, butane, pentane and C<sub>6</sub>+ (hexane and higher hydrocarbons)

#### OFFICIAL METHODS MET:

ASTM 1945, ASTM 1946, GPA 2177;


#### KEY FEATURES & BENEFITS:

Gas analysis is performed using this complete system that allows to determinate H<sub>2</sub> using the channel equipped with a  $\mu$ TCD, a micro-packed column and a micro-packed column, an 8-ports valve and a solenoid valve for cleaning the system by other compounds;


Permanent gases and hydrocarbons analysis is performed using the channel equipped with the  $\mu$ TCD, micropacked columns, a 10-ports valve and a 6-ports valve to bypass the Molsieve column.

The gas sample is introduced separately in the single analytical channels by the 8-ports and 10-ports valves equipped with sampling loops to 1 ml.


All the analytical parameters can be controlled by the DDS Clarity™ Chromatography Station Software;



Chromatogram of the gas standard mixture on Channel 1  
The peak corresponds to He+H<sub>2</sub>



Chromatogram of the gas standard mixture on Channel 2.  
The peaks correspond to: CO<sub>2</sub>, C<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, C<sub>3</sub>, iso-C<sub>4</sub>, n-C<sub>4</sub>, iso-C<sub>5</sub>, n-C<sub>5</sub>



Chromatogram of the gas standard mixture on Channel 2.  
The peaks correspond to: CO<sub>2</sub>, C<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>

**EXPERIMENTAL PARAMETERS FOR PERMANENT GAS, CO<sub>2</sub>, HYDROCARBONS AND SULPHUR COMPOUNDS**

Master GC Analyzer	
Columns	Channel 1 <b>Hayesep Q</b> (1m, 1/16", 80/100 mesh) <b>Molsieve 5A</b> (2m, 1/16", 80/100mesh)
	Channel 2 <b>Hayesep T</b> (1m, 1/16", 80/100mesh) <b>Hayesep Q</b> (2m, 1/16", 80/100mesh) <b>Molsieve 13X</b> (2m, 1/16", 80/100mesh)
GC Oven	
Temperature	Temp: 32°C, Time: 23 min, Rate: 10°C/min Temp: 32°C, Time: 13 min
Injector A: PK	
Temperature	OFF
Carrier Gas	Nitrogen
Column Type	Capillary - 200m, 0.53mm
Flow	Set: 14 ml/min - Measured at Det C: 5.2 ml/min
Injector B: PK	
Temperature	OFF
Carrier Gas	Helium
Column Type	Capillary - 200m, 0.53mm
Flow	Set: 15 ml/min - Measured at Det C: 2.8 ml/min
Inj C: PK	
Temperature	OFF
Carrier Gas	Helium
Column Type	Capillary - 200m, 0.53
Flow	Set: 13 ml/min - Measured at Det C: 2.7 ml/min
Det A : μTCD Channel 1	
Temperature	160°C
Main Filament Temperature	190°C
Filament Safety	Injector A
Min. Half-Peak Width	0.6 s
Digital Acquisition Rate	25 Hz
Signal Zeroing	10 mV
Det A : μTCD Channel 1	
Temperature	160°C
Main Filament Temperature	190°C
Filament Safety	Injector A
Min Half-Peak Width	0.6 s
Digital Acquisition Rate	25 Hz
Signal Zeroing	10 mV
Det C : μTCD Channel 1	
Temperature	160°C
Main Filament Temperature	190°C
Filament Safety	Injector C
Min. Half-Peak Width	0.6 s
Digital Acquisition Rate	25 Hz
Signal Zeroing	10 mV
Auxiliary Pressure	
Aux Press	2.2 bar
Auxiliary Temperature	
Aux Temp (°C) - AUX Oven Temp	100°C



## REFINING GAS ANALYSIS

По вопросам продажи и поддержки обращайтесь:

Архангельск (8182)63-90-72  
 Астана (7172)727-132  
 Астрахань (8512)99-46-04  
 Барнаул (3852)73-04-60  
 Белгород (4722)40-23-64  
 Брянск (4832)59-03-52  
 Владивосток (423)249-28-31  
 Волгоград (844)278-03-48  
 Вологда (8172)26-41-59  
 Воронеж (473)204-51-73  
 Екатеринбург (343)384-55-89  
 Иваново (4932)77-34-06  
 Ижевск (3412)26-03-58  
 Казань (843)206-01-48

Калининград (4012)72-03-81  
 Калуга (4842)92-23-67  
 Кемерово (3842)65-04-62  
 Киров (8332)68-02-04  
 Краснодар (861)203-40-90  
 Красноярск (391)204-63-61  
 Курск (4712)77-13-04  
 Липецк (4742)52-20-81  
 Магнитогорск (3519)55-03-13  
 Москва (495)268-04-70  
 Мурманск (8152)59-64-93  
 Набережные Челны (8552)20-53-41  
 Нижний Новгород (831)429-08-12  
 Новокузнецк (3843)20-46-81

Новосибирск (383)227-86-73  
 Омск (3812)21-46-40  
 Орел (4862)44-53-42  
 Оренбург (3532)37-68-04  
 Пенза (8412)22-31-16  
 Пермь (342)205-81-47  
 Ростов-на-Дону (863)308-18-15  
 Рязань (4912)46-61-64  
 Самара (846)206-03-16  
 Санкт-Петербург (812)309-46-40  
 Саратов (845)249-38-78  
 Севастополь (8692)22-31-93  
 Симферополь (3652)67-13-56  
 Смоленск (4812)29-41-54

Сочи (862)225-72-31  
 Ставрополь (8652)20-65-13  
 Сургут (3462)77-98-35  
 Тверь (4822)63-31-35  
 Томск (3822)98-41-53  
 Тула (4872)74-02-29  
 Тюмень (3452)66-21-18  
 Ульяновск (8422)24-23-59  
 Уфа (347)229-48-12  
 Хабаровск (4212)92-98-04  
 Челябинск (351)202-03-61  
 Череповец (8202)49-02-64  
 Ярославль (4852)69-52-93